

## TRNAU1AP-tRNA selenocysteine 1 associated protein 1 Gene

Reference: PTXBC039879

|                |                                                                                                                                                                                                             |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                         | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                        | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | PRO1902; SECP43; TRSPAP1; tRNA selenocysteine 1-associated protein 1; tRNA selenocysteine associated protein (SECP43); tRNA selenocysteine-associated protein 1; tRNA selenocysteine 1 associated protein 1 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | tRNA selenocysteine 1 associated protein 1                                                        |
| <b>GENE ID</b>          | 54952                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC039879                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | TRNAU1AP                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgctgtatgaattcttctcaaaagtctacccctctgtcggggaggcaaggtgttttgaccagacaggcgtgtctaagggttatgtttgtgaaatcacagatgaactggaacagaagcgagccctgacggagtccagggagcagtgaggactgggtctaaagcctgtcggctgag  
cgtggcaatccctaaagcgagccgtgttaaagccagtggaatagtcagatgtacagttatagctacaaccagtattatcagcagttaccagaactactatgctcagtggggctatgaccagaacacaggcagctacagctacagttacccccagttatgctataccagagcaccatgcaga  
catatgaagaagttggagatgatgcatfaggagcccatgccacagctggatgtgactgaggccaacaaggagttcatggaacagagtgaggagctgtatgacgctctgactgactggcagccctggacacagtgctctcagagatccctgcatgatgtag