



Reference: PTXBC039904

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AHO3; BDMR; HA6116; HD4; HDAC-4; HDAC-A; HDACA; histone deacetylase 4; histone deacetylase A		

For research use only.

GENE DESCRIPTION	histone deacetylase 4
GENE ID	9759
ORIGIN SPECIES	Human
ACCESSIONS	BC039904
VECTOR	pDONR223
NCBI SYMBOL	HDAC4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]