

C1orf110-chromosome 1 open reading frame 110 Gene

Reference: PTXBC040018

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf110; coiled-coil domain-containing protein 190; coiled-coil domain containing 190		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 110
GENE ID	339512
ORIGIN SPECIES	Human
ACCESSIONS	BC040018
VECTOR	pDONR223
NCBI SYMBOL	C1orf110
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtcaggggacaaattgtataagcattttgattggagagggaatgccagcaggctgaagccagactggacaaagactgcagagactaaaggttattgctctaccatgtgaaattgctgacctgggagcagaggcagctccaaaagaactgcagaggttcagcagcagaa
accatgaagaaaaagttctctcttatttgggaatggatttcagaagagaccagaagatgttctctgttctcaccacaggggaagcagaagcacagagcccccacaggcctaagaaaatgagagcattggcaaccctgtatggccaagacacatgcaaaagcaagtcagggtcctcct
tcacatgatgctggcctcaagaccccatgaagagcaaaaagcagccactctcctcaaaataacagaactgctgcttcataaaagagcaaccacaagcccaagagaaagattctgtgaatccatcaaggacgtagacccagcaagggcatctctgttccatgccaaatcaagaggtt
ccaccaacacatagaacaaggtccttagtccagccagctagtattggaatggcatgtgctga