

CRCP-CGRP receptor component Gene

Reference: PTXBC040107

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGRP-RCP; CGRPRCP; RCP; RCP9; DNA-directed RNA polymerase III subunit RPC9; CGRP-receptor component protein; RNA polymerase III subunit C9; calcitonin gene-related peptide-receptor component protein; CGRP receptor component		

For research use only.

GENE DESCRIPTION	CGRP receptor component
GENE ID	27297
ORIGIN SPECIES	Human
ACCESSIONS	BC040107
VECTOR	pDONR223
NCBI SYMBOL	CRCP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaagtgaaagatgccaaattctgccttctcagtaactacgagacgttaaaatacatatcaaaaacacatgcaggcaccagagtcctgaaattgtcagagaattctcacagcattgaaaagccacaagttgaccaaagctgagaagctccagctgctgaaccaccggcctgtgactgctgtggagatccagctgatggtggaagagagtgaaagcggctcagcggaggagcagattgaagctcttctccacacctcaccagcattctgcctgcagagccagaggctgagcagaagaatacaaacagcaatgtggcaatggacgaagaggacccagcatag