

MED22-mediator complex subunit 22 Gene

Reference: PTXBC040111

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MED24; SURF5; surf-5; mediator of RNA polymerase II transcription subunit 22; surfet locus protein 5; mediator complex subunit 22		

For research use only.

GENE DESCRIPTION	mediator complex subunit 22
GENE ID	6837
ORIGIN SPECIES	Human
ACCESSIONS	BC040111
VECTOR	pDONR223
NCBI SYMBOL	MED22
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccagcagagagccctgccccagagcaaggagacgtgctgctgagtcctacaacaagcggtgaaggacgacattaagtcctcatggacaactcaccgagatcatcaagaccgccaagattgaggacgagacgcaggtgtcacgggccactcagggtgaacaggacaatta
cgagatgcatgtcgagccgccaacatgctccgagccgagtcctgatgaagctggtgtccgacctcaagcagttcctgacctcaatgacttccctcgtgaacgagccattgaccagcgcaaccagcagctgctgcacactgcaggaggagtcgaccggaagctcatcag
ctcgagagcagatctccattgacctctacgagctggaggaggagtattactgtccaggtataaataag