

TMEM9B-TMEM9 domain family, member B Gene

Reference: PTXBC040124

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C11orf15; transmembrane protein 9B; TMEM9 domain family member B		
For research use only.			

GENE DESCRIPTION	TMEM9 domain family, member B
GENE ID	56674
ORIGIN SPECIES	Human
ACCESSIONS	BC040124
VECTOR	pDONR223
NCBI SYMBOL	TMEM9B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgacccctgtgggaggcctctcggccttgcctcctgtcagcctgtcgtcctggcgttccgtgctgctggcgagctgtcagacggcgccaagaatttcgaggatgtcagatgtaaatgtatctgccctccctataaagaaaattctgggcattttataaagaacatatctc
agaaagattgtgattgcctcatgttggagccatgcctgtgcggggcctgatgtagaagcatactgtctacgtgtgaatgcaaatatgaagaaagaagctctgtcacaatcaaggtaccattataattatctctccatttggcccttctactctgtacatggtatatcttactctggtgag
cccatactgaagaggcgctcttggacatgcacagttgatacagagtgatgatgataatggggatcaccagcctttgcaaatgcacacgatgtgctagcccgtctcccgagtcgagccaacgtgctgaacaaggtagaatatgcacagcagcgtggaagcttcaagccaagagcagc
gaaagtctgtcttgaccggcatgtgtcctcagctaa