

DOLPP1-dolichyl pyrophosphate phosphatase 1 Gene

Reference: PTXBC040152

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LSFR2; dolichyldiphosphatase 1; dolichyl pyrophosphate phosphatase 1; linked to Surfeit genes in Fugu rubripes 2		

For research use only.

GENE DESCRIPTION	dolichyl pyrophosphate phosphatase 1
GENE ID	57171
ORIGIN SPECIES	Human
ACCESSIONS	BC040152
VECTOR	pDONR223
NCBI SYMBOL	DOLPP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctccagccattccagtttatgtgttctctccgtctattctctctttctctgtatttaagaatgcaccaacaacaacgccaggttctctggacttgcgtgtgaggcagctgtctctcctgggactctcgcgtgtggccttctagtctctacagcaggcctgtctccgagttcttcta
ccgagacacaagcctattcccaacgtactctgtttgagtacacggtaacccggcagaaagccaggaacagacaacgcaagctggggacgaaactgcagtga