



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C13orf33-chromosome 13 open reading frame 33 Gene

Reference: PTXBC040165

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

C13orf33; AWMS3; MEDA-4; MEDA4; hAWMS3; mesenteric estrogen-dependent adipogenesis protein; activated in W/Wv mouse stomach 3 homolog; mesenteric estrogen-dependent adipose 4; mesenteric estrogen dependent adipogenesis

For research use only.

GENE DESCRIPTION	chromosome 13 open reading frame 33
GENE ID	84935
ORIGIN SPECIES	Human
ACCESSIONS	BC040165
VECTOR	pDONR223
NCBI SYMBOL	C13orf33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctggggggcgccctgcgagcgtgtgtgccagggcgagcctgacctcaatctcgtctggggagcttcgcagcctgtggagcttgcgactggcgagctggccctctgcgcctgttcagctcgtctgcgcctgcagcccggtgccttcacgtgagcggcgaccacgtctgtggfgccgg
gcccggggagcggcgccggcgccgggggggcttcaacgtcttcggtgacggcctctgctggcctcgacggcgagctctaccgctcagcagctacatcaagagtgatgtggaactgaccaactactgtgattataaagactacagggaaactatattgagcaaaccaattgtgtcttatt
aatgtacagacacaaaagacacctcaaaagaagggacgtacgctgtttctgtaaacggcgaccgccccaagataaagaagacagatagagaaggatggcagcatggatctctcctcagtgattggagaagtattaccggcttcagtttgattttcaagagcagtgagaagtattctcccc
aggaataagaatgttgtaaggaaaatttaagcttgcataatgaattcaacgtctgatctcaattattgatttcttatttggtttggcctcgaattcctgtgtgaaaggtaaatgaaaaggtaattctgaatttgcaggtacaagctccagaagaagagacgattcaatttctggaataaagtgtg
ccttcaattcaagaaagcagcaattcttcgcagaaattcagtttaacttcacagaggtttcaalagatgaatttttgcctcaactctcagccagcatcagctatctgcagagcggcttttggcagataccatttgcgagatcaaccattcaagtttctgaattctgaatctgaatctgaatctgaatctga

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1