



Reference: PTXBC040168

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ERG, ETS transcription factor; transcriptional regulator ERG (transforming protein ERG); erythroblast transformation-specific transcription factor ERG variant 10; TMPRSS2/ERG fusion; transcriptional regulator ERG; erg-3; p55; ets-related; v-ets avian eryt		

For research use only.

GENE DESCRIPTION	v-ets erythroblastosis virus E26 oncogene homolog (avian)
GENE ID	2078
ORIGIN SPECIES	Human
ACCESSIONS	BC040168
VECTOR	pDONR223
NCBI SYMBOL	ERG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]