



Reference: PTXBC040239

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	pyruvate dehydrogenase kinase, isozyme 4; [Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 4, mitochondrial; [Pyruvate dehydrogenase [lipoamide]] kinase isozyme 4, mitochondrial; pyruvate dehydrogenase kinase, isoenzyme 4; pyruvate dehydrogen		

GENE DESCRIPTION	pyruvate dehydrogenase kinase, isozyme 4
GENE ID	5166
ORIGIN SPECIES	Human
ACCESSIONS	BC040239
VECTOR	pDONR223
NCBI SYMBOL	PDK4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1