

MRPL42-mitochondrial ribosomal protein L42 Gene

Reference: PTXBC040240

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC204; L31MT; L42MT; MRP-L31; MRP-L42; MRP-S32; MRPL31; MRPS32; PTD007; RPML31; S32MT; 39S ribosomal protein L42, mitochondrial; 28S ribosomal protein S32, mitochondrial; 39S ribosomal protein L31, mitochondrial; mitochondrial ribosomal protein L42		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L42
GENE ID	28977
ORIGIN SPECIES	Human
ACCESSIONS	BC040240
VECTOR	pDONR223
NCBI SYMBOL	MRPL42
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtagctgcagtaaaatgggtgatgtcaaaagagaactatcttgaacatttattccagtcctcaaaatggagctttatattgtgtttgtcataaaatctacgtattctctctaccagatgactataattgcaacgtagagcttgcctgacttctgatggcaggacaatagtagctaccaccttctgtggacattccatgaacacacaaacatccctcgccagatcctgtgcataataatgaagaacacatgatcaagtgcctgaaaaccagattggaagaaaaagtgaacaccttgaggagggacatgatagaacaacttagcaaaatgttctttactactaagcaccgttggatcctcattggacggtatcacagatgtcgaagaatctgaatcctccaaaagacagatga