

## SFRS11-splicing factor, arginine/serine-rich 11 Gene

Reference: PTXBC040436

|                |                                                                                                                                                                                                                                                 |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | SFRS11; NET2; dJ677H15.2; p54; serine/arginine-rich splicing factor 11; SR splicing factor 11; arginine-rich 54 kDa nuclear protein; splicing factor p54; splicing factor, arginine/serine-rich 11; serine and arginine rich splicing factor 11 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | splicing factor, arginine/serine-rich 11                                                         |
| <b>GENE ID</b>          | 9295                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC040436                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | SFRS11                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgagcaacactaccgtctccccagcactgcaggctccggggcccccagggcgccggcggtggcgagggtgggtggcgccggaggcgccggcaccgaggttaaccaggtgactaatgtctccccgagcgctagctctgagcagatgcggactctctcggttcttaggcaagatc  
gacgaactgcgcctcttcccgccggatgattcgcttccagctctcatctcgtgtctgtttgtaagtccatgatccagactcagcagttgtggcacagcatctgacaaacactgtattcgttgacagagctttgatagtcgtaccatgacagaaggagttattctgatgaagctaaagctttg  
tctctgttgccaccagctaatgcagtgccaggcttctgctgtgtggactcctgcctactcctaaccactaccagattggcgctgtccactggctgtttggggcgctcctactcttgatcctgccccttgcgtcacttgggcttctggagcaaaccttgaactctcagctcttgcgcagat  
cagttgctgaagcttatgagtactgttgatcccaagttgaatcatgtactgtgctgtctgtttcaccaagctcgaatcggataccctctagtaaaagaataagaggaagctatgaaaagagtagcagaagcagtccttaatttctgctgctatagaaccagataaagaagaagaaaaagaa  
ggcattcaagatcaagatcacgttctaggaaggaggagactccctcatctctagacacaggcggtcaagaagcagatcgagacggcggtcacattctaaagtctaggaagtcggcgacgatccaaagcccaaggcgagagaagatctattccagagaagaggtagaagggtcaaggaa  
gcacatcaaaaacaagagacaaaaaagaagaagacaagaagaagaacgttctaaacaccaccaaaaagtacagcacagccagacgttctagaagtgcaagcagagagagacgacgacgaagaagcaggagtgccacaagatctctaaaaagcctcggtctctaaaaagaaa  
attgtcccgctcaccatcccttaggagacataaaaaaggagaagaagaagataaagacaaaagaaagtagtgggagtagaagagaacgatcaacaagcaagaagaagaagtagaagataaaggaaagacggggaagaaaaatcagagagtgataaagatgtaaaagttacacgg  
gattatgatgaaggagaacagggtgatgacagtgaagaagagaaaaagaagaagaacaaatagaacaggttcccttaaaacaagaagatgtctgtgaaaagggaactggtgattcactaagagaatcacaagtgatggggtgatcatatgaagaagacatgatatga  
gtgactga