



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

AQP11-aquaporin 11 Gene

Reference: PTXBC040443

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

AQPX1; aquaporin-11; AQP-11; aquaporin 11

For research use only.

GENE DESCRIPTION	aquaporin 11
GENE ID	282679
ORIGIN SPECIES	Human
ACCESSIONS	BC040443
VECTOR	pDONR223
NCBI SYMBOL	AQP11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcgccgctgctggggctccggtccgagctgcaggacacctgcacctgcctggagctgatctgctgggtgctgctatcggggctggcccgctagtcgcccgagcagctgcacagccggctgggccacgctctgctctggagttctagccacctccagctctgctgctc
accacagagctgcacaactctgtagcgaacagcaccgccgaccaccactggagctgagctctgctactcttctccttgccttggctgatctgcttggggcagctccagcaaccgctgcggcgtgatgatcagatgatctggggggcatgtcccccagacgggctgggga
ggctattggctcagctgcttagtgcctctgtagcagctactgcacaagcgccttggagcttggctgctgaccagatacactgcagcagagaggctctgccttgcagaataccatccagctgacttgcctcaagcggctcatcacagaggccgtctgctctcttctctccacagcgtct
gctgcaacttcagcgaagtcgaaccaagctctgtaaccacctgctggctgcactatccacttttggctctatgcaggagggaagcttaacagagctgtattataccgcttggcctgctcatcttcatctttgatgaagcattccctcagtttttatgatactgctgctctcttta