

IL23R-interleukin 23 receptor Gene

Reference: PTXBC040720

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	interleukin-23 receptor; IL-23 receptor		
For research use only.			

GENE DESCRIPTION	interleukin 23 receptor
GENE ID	149233
ORIGIN SPECIES	Human
ACCESSIONS	BC040720
VECTOR	pDONR223
NCBI SYMBOL	IL23R
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaaaacagcaatgtgtgaaaaatgctacagaaaaatggaacttatgaataaattccagtgcagcaggtcctatatgttgatcccatgattacagagataaaagaaatcttcatccagaacacagcctacagactacaagaaggagaatacaggaccctggagacaagagactac
ccgcaaaactcgtattcgacaatactacagttgtatataatctctgatctcaacactggatataaaccccaatttcaaatgttgcctgagggagccatctcagcaataataatgaaattacttcttaacacttaaacaccagttgattccttagactcaggaaataatccaggttacaagaac
atcctaatttgccttttctgttcaagtgtgaattcactaagcaacacaataattcttggagaattagccctcatattaaatcaaggagaatgcagttctctgacatacaaaactcagtagaggaggaaaccacatgcttttggaaatgattcaccagtgaaactattccagaacagaccctgct
tcctgatgaattgtctcctgtttggggatcgtgaatgaggagttgccatctattaatactattttccacaaataatttggaaagccacttcaataggatttcactcttggaaaagtag