

## KLK8-kallikrein-related peptidase 8 Gene

Reference: PTXBC040887

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HNP; NRPN; PRSS19; TADG14; kallikrein-8; ovasin; serine protease 19; serine protease TADG-14; tumor-associated differentially expressed gene 14 protein; kallikrein related peptidase 8		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	kallikrein-related peptidase 8
<b>GENE ID</b>	11202
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC040887
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	KLK8
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

```
atgggacgccccgacctgtgcggccaagcgtggatgttctgtcttggggggagcctgggcaggacactccagggcacaggaggacaagtgctgggggtcatgagtgccaaccccttcgagccttggcaggcggcctgttccaggccagcaactctgtggc
gggtgctctgtaggtggcaactgggtccttacagctgccactgtaaaaaacgaaatacacagtacgcctgggagaccacagcctacagaataaagatggcccagagcaagaataacctgtgttcagtcctcccaacccctgtacaacagcagcgatgtggaggaccacaacct
gatctgatgcttctcaactgctgaccaggcatccctgggtccaagtgaagccatcagcctggcagatcattgcaccagcctggccagaagtgcaccgtctcaggctggggcactgtcaccagtccccgagagaatttctgacactcacaactgtcagaagttaaaatcttccc
cagaagaagtgtgaggatgcttaccggggcagatcacagatgtcatgtctgtcaggcagcgaaggggctgacacgtgccaggcgtattctggagcccccctggtgtgtgatgtgcactccaggcatcacatctgggtcagaccctgtgggaggtccgacaacct
ggcgtctataccaacctgtccgctacctggactggatcaagaagatcatagcagcaaggcctga
```