

CWC15-CWC15 spliceosome-associated protein homolog (*S. cerevisiae*) Gene

Reference: PTXBC040946

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CWC15 spliceosome associated protein homolog; CWC15 spliceosome-associated protein; spliceosome-associated protein CWC15 homolog; AD002; C11orf5; Cwf15; HSPC148; ORF5		

For research use only.

GENE DESCRIPTION	CWC15 spliceosome-associated protein homolog (<i>S. cerevisiae</i>)
GENE ID	51503
ORIGIN SPECIES	Human
ACCESSIONS	BC040946
VECTOR	pDONR223
NCBI SYMBOL	CWC15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacaacagcagccaggccaaccttgaacctgccagaggtggaaggggaaaaaggagaaggtgattgagccaacttcaaagcagattcaagcagagacctaccctctatacaaaataacagacagactactcaggatgccctgaagaggttcgtaacctgactcagg
agagagttggaagaagagagagagctgctgcaagagagaaaaatagggatcgtccaacccgagaacatacaacctctcttcagtggtcaaaaaagccacgggttagaccagattctgccgccaaccttgatgcagatgacctctaacagatgaggaagatgaagatttgaagaaga
aagtgatgatgatactgcagctctcttcgagaactggaaaaaattaaaaaagaagagctgaagagcaggccagggaaggaacaagaacaaaaagctgaagaagagagattcgtatggaaaacattctgagcggaaacctctcttaatctcactggccatccagcctcaggc
caacttcaaggttaaaagaaggtgggatgacgttcttcaagaactgtgcaaaaggtgtagatgaccagaagaaagacaaaagattgtaaatgacacactgcgactgaattcacaaaagttcatggagaatatattaaatag