



Reference: PTXBC041028

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CARKD; LP3298; ATP-dependent (S)-NAD(P)H-hydrate dehydratase; ATP-dependent NAD(P)H-hydrate dehydratase; ATP-dependent NAD(P)HX dehydratase; carbohydrate kinase domain containing; carbohydrate kinase domain-containing protein; NAD(P)HX dehydratase		

For research use only.

GENE DESCRIPTION	carbohydrate kinase domain containing
GENE ID	55739
ORIGIN SPECIES	Human
ACCESSIONS	BC041028
VECTOR	pDONR223
NCBI SYMBOL	CARKD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr