

CCDC103-coiled-coil domain containing 103 Gene

Reference: PTXBC041060

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CILD17; PR46b; SMH; coiled-coil domain-containing protein 103; coiled-coil domain containing 103		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 103
GENE ID	388389
ORIGIN SPECIES	Human
ACCESSIONS	BC041060
VECTOR	pDONR223
NCBI SYMBOL	CCDC103
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaggaatgacatcatcaactcaaggcttggagaaagagctgaggctgactcactgctgatgagaagtacaacgggagaatgctgccaagttacggcagtggaacagagggtggcttctatgaggagttcaggggtattgcttgcacatctgaagccactggagc
ggaaggataagatgggaggaagagaactgtgcccgtggaactgtcacactattcagggaaggaccttcaggatgtggccactgaaatcctccggagaaagccccctccagcccgagacgtctgctgacttctatctgattggcgacgacactgccaagtggccagagcgtac
caggctctactgcagcttgggggtccaaggctcggctgctcttcagacagatgtgggattggacttctggggagctgctgtggcactggctgacgtggggccggctgaccgggcagcgggtgctggggatcctatgcagcctggcgagcactgggcgttcacctgaaccta
agcctgctgagccgggcagagagagagagctgcaagggtgtttcagaagctgcaagccatgggcaacccagatccgtgaaggaggggctcagctgggaggagcaggggtctggaggagcagctgtgtgggtccaggaggagagaggctcctcaggagctgctggagct
gtaccaggttgactga