

ZNF695-zinc finger protein 695 Gene

Reference: PTXBC041082

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SBZF3; zinc finger protein 695; zinc finger protein SBZF3		
For research use only.			

GENE DESCRIPTION	zinc finger protein 695
GENE ID	57116
ORIGIN SPECIES	Human
ACCESSIONS	BC041082
VECTOR	pDONR223
NCBI SYMBOL	ZNF695
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggactattggcattcaggatgtggctctagaattctccagaggagtggaatgcctggaccagctcagcggagttgtataggatgtgttagagaactacagaaacctgatcccttggtgaggatagctcaatgcattcctatttcacagtcttgctatgtaagccag
aactgatcatctgtctggaggcaaggaaagagccctggaactgaaacacagagaagacagccaaacactcagcttctctcttactgaagacattttccagagcagggcctgcaagtttcattccaaaagtgtgctgagaagatgaaagatgttctctgagaattacgcta
aggaatgactgggaaattccatgtgaagatgtgcttgccttcccttgccttctgcatgattctaagtttctgaggcctccagaaagcagaagcatgtaaagccacagaaaccggtccagtcgtcgaccatcgagctgcttga