



Reference: PTXBC041091

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WHITE2; ATP-binding cassette sub-family G member 4; ATP-binding cassette, sub-family G (WHITE), member 4; ATP binding cassette subfamily G member 4		

For research use only.

GENE DESCRIPTION	ATP-binding cassette, sub-family G (WHITE), member 4
GENE ID	64137
ORIGIN SPECIES	Human
ACCESSIONS	BC041091
VECTOR	pDONR223
NCBI SYMBOL	ABCG4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]