

TAF5L-TAF5-like RNA polymerase II, p300/CBP-associated factor (PCAF)-associated factor, 65kDa Gene

Reference: PTXBC041094

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PAF65B; TAF5-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 5L; PAF65-beta; PCAF-associated factor 65 beta; TAF5-like RNA polymerase II, p300/CBP-associated factor (PCAF)-associated factor, 65kDa; TATA-box binding prote		

For research use only.

GENE DESCRIPTION	TAF5-like RNA polymerase II, p300/CBP-associated factor (PCAF)-associated factor, 65kDa
GENE ID	27097
ORIGIN SPECIES	Human
ACCESSIONS	BC041094
VECTOR	pDONR223
NCBI SYMBOL	TAF5L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaacgagtgctaccgagcagattcagatggcagtgctctgctacccaacgcgcgagtgactgagtgctccctgaagcaaggactgctgctacagactgctgaagagatggcgccaatcaacagtgcaatcagaatctggtgtgccacatagtgctgca
gccccctggccaggcagaacccagcaatgaaglacagttggagcactgcggaattttctcactgattctgattccagcatagcccaagagtgatgctctctctatctctctttgtctacccatctcaacctgggtccaaaacagtcggaagagcacagtggaaagtgttacagccgct
tccatggaatgtttctgcagaatgctagccagaaggatgctattgagcagctacagaccactcaaacatccaggacatctctatctaaactcaagcttcgagcattctagataacaagtagctggtgctcctccaaagagacagctacaactaccttaccgctacccaaagtgaacaat
actgccccgtgcaaaagtcctacattacatattctgacgtgcagcctgccaagagaacagactatcagctgtatgccagtgagcagctctcccgagtgagaacaacgggttgagccccctgacatgccagccctattctgcagaacgaggtgcccctagaggtcttacaggagag
cattaagcgcgtcaaggatgggctcctccctcactaccatctgctctatgctctctataacacagagcagctgttgaaacactgcagaatctccccgatagcaagctgctgctgctgggttgacaactcctgtataaaactttggagttacgatccaagaagttaaatcagagcccca
ccaagtagacgtgtcccgcatcattggcctgtgatattctggaggaggaggtatga