

MRPL43-mitochondrial ribosomal protein L43 Gene

Reference: PTXBC041165

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L43mt; MRP-L43; bMRP36a; 39S ribosomal protein L43, mitochondrial; mitochondrial ribosomal protein bMRP36a; mitochondrial ribosomal protein L43		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L43
GENE ID	84545
ORIGIN SPECIES	Human
ACCESSIONS	BC041165
VECTOR	pDONR223
NCBI SYMBOL	MRPL43
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacggcgccgggactccgagccgtcttggccagcgttctccacaacggactgggtcgtatgtgcagcagctgcagcgtctgagctcagcgtcagccgcgacggcgccctcgtctcgcggcgccaggaggtctgtgagcgggaggtgatcgaactcgcggcaggaatcca
ggggtcgtatataatgtaaacctcgcgtcgtcgtgcccagagtagtggccgaataaccttaacgggctgtgcgcgaggagagatccactgcaagtcggtcagggagatctgacgctggtgcagaagctggccgaccagtcgggcttgacgtgatccgacatccgaagccctt
ccacaccgacaaccctagcatccaggccagtggtcacccttcaccaacaagccgaccacgttccgctgggtacgccccgagaggttcaggatcctgccccagccaggacactggcctgagactgtctgcagttgcacgcagatcctctgccccgtggccgaccaccag
acctccccacagtgatcctatctcatcctcattgacctctgtccagccccctatgtctccgagtttctgctcccgattgtccctgcactgaccactgtgtctcagcgtga