



Reference: PTXBC041167

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TBC1 domain family member 14		
For research use only.			

GENE DESCRIPTION	TBC1 domain family, member 14
GENE ID	57533
ORIGIN SPECIES	Human
ACCESSIONS	BC041167
VECTOR	pDONR223
NCBI SYMBOL	TBC1D14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr