

C22orf25-chromosome 22 open reading frame 25 Gene

Reference: PTXBC041339

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C22orf25; MECRCN; transport and Golgi organization protein 2 homolog; transport and golgi organization 2 homolog		

For research use only.

GENE DESCRIPTION	chromosome 22 open reading frame 25
GENE ID	128989
ORIGIN SPECIES	Human
ACCESSIONS	BC041339
VECTOR	pDONR223
NCBI SYMBOL	C22orf25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcatcattcttaagttgatcctgcctgtttccaaaacgcgtacagctcatcttgccagccacaggatgaattctacagccgacctccaagttagctgacttctgggggaacaacaacagatcctcagtgaggctggacatggaggaaggcaaggaggacacatg
gctgggcatcagcacacgtggcaagctggcagcactcaccaactacctgcagccgcagctggactggcagggccggagggcgaggtgaactgtcaccactttctgaccactgacgtggacagctgtgtctacctaagaaggctctctatggagggccatctgtacaatggctcaacct
catagcagccaacctgagcacagcaaggagacgtcatttgcatactatgggaaccgaggggagcctgacatcgttttgacgccaggcacctacgggctgagcaacgcgtctgagactccctggaggaagctgtgctttgggaagcagctcttctggaggtgtggaacgga
gccaggcgctgcccaaggatgtgtctatcgcagcctcctggatgtgtctcaacaataaaggcgagctgccagaccggccatcgaggaccagggtggggagtacgtgcagccatgctgagcaagtlacgaggctgtgtgctgtgcgtgcctggctacggcaccagaaccaac
actatcatcctgtagatgcggagcggccacgtgacctcactgagcgttagcatgatggacaaggacctctcccactggagaccagaacctatgattcacactgcagagctaa