

## FRMD3-FERM domain containing 3 Gene

Reference: PTXBC041376

|                |                                                                                                                                                                                    |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                               | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | 4.1O; EPB41L4O; EPB41LO; P410; FERM domain-containing protein 3; band 4.1-like protein 4; band 4.1-like protein 4O; ovary type protein 4.1; protein 4.1O; FERM domain containing 3 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | FERM domain containing 3                                                                         |
| <b>GENE ID</b>          | 257019                                                                                           |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC041376                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | FRMD3                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgatccatttcggagctccagctcaaatcgctcagccaggagatgagatgcacatccggctgctggagcactcggagatctcctgccatccagagggaacaaagggcagtttctcattgaccacatctgcaactactacagcctgctggagaaggactattggcattcgeta  
tfgggaccagagaagcaaaaggcactggcttgaacctaacaagtccatctcaagcaaatgaaaactcatcccatcacacatgtgcttttagagtgaattctaccacatgaaccttgaagattaaagaagagctcacaagatacctttatcttcagattaaaggacattttcatgg  
ccactctgtgctcctttctgatgctcctacctgggtgctgtattgttcaagctgagctgtggtgattacgacccaatgagcatcctgagaaattacatcagtgagtttgagatttccccaaagcagtcacagaagctggaaagaaaaatagtggaattcataaaatgaactcagggcgag  
agcccacagttgctgaatttaactgtcctgaaagctcacacttggaaacctacgggtggatcctcaccatgcaaggattcaacaggcacaacaacatttttaggattcacagctcaggtcttgggtctttcagggaataagagaatccattgataaaatggccagatgtctcaaa  
ttgaagtttgaagggaagacattttatgtgattggcaccagaaggagaaaaagccatgttggcattcaatacttcaacaccagctgcttgcacacatcttggaaagtggagtggaagaaaccaggccctttataagtagtcaaaatccagtcagatcaagactgtatcaagcagcaagatattt  
tttaagggaagtagatttcgatatagtggaaagtggccaaagaggtggtggagccagttccaagatccagaggagcctcctgaggtgcacagagccaactacagagccgagctccactccttgaacaaacagctcatcattaaatggaacccctgcagccctgcttcttc  
ccccagcagcagaagaagaactcctctgggtgaggtgttccattgcttaagaggagaacatttctgctccttgcctcagctccagtgaaagcagccgggagatgaagatccccctagtgaaaggaagataaaataaagaagaacatttaaccatctgaactagtgt  
acaacccaagtccagcctgctccccccccctgtggatgacgatgagattgacatgctcttgcctctctagcttgccttggaaagagaagacacagattcattgagatctggaagcagatgaaaacgaccttttgcctgctgaagaagaggagctgaaggagctcgcctgctt  
gtcgtggagctacgacattctgactggccatattcgggtgaacccactgtcaagagttttccaggctccttgggtggcctgggactgctcttcttatttccccctgctcctcctcttggagtcaggtattgatctctcttctatgcgaaatccgcagacaccagagtttgagcagttt  
cactatgaatactactgtccctcaaggagtggtggtgggaaagtcacatcctctacatgctgggttgcctatga