



Reference: PTXBC041407

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GLUR3; GPRC1C; MGLUR3; mGlu3; metabotropic glutamate receptor 3; glutamate receptor, metabotropic 3; glutamate metabotropic receptor 3		

For research use only.

GENE DESCRIPTION	glutamate receptor, metabotropic 3
GENE ID	2913
ORIGIN SPECIES	Human
ACCESSIONS	BC041407
VECTOR	pDONR223
NCBI SYMBOL	GRM3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr