



Reference: PTXBC041560

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MOAT-B; MOATB; multidrug resistance-associated protein 4; MRP/cMOAT-related ABC transporter; bA464I2.1 (ATP-binding cassette, sub-family C (CFTR/MRP), member 4); canalicular multispecific organic anion transporter (ABC superfamily); multi-specific organic		

For research use only.

GENE DESCRIPTION	ATP-binding cassette, sub-family C (CFTR/MRP), member 4
GENE ID	10257
ORIGIN SPECIES	Human
ACCESSIONS	BC041560
VECTOR	pDONR223
NCBI SYMBOL	ABCC4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]