

FARP1-FERM, RhoGEF (ARHGEF) and pleckstrin domain protein 1 (chondrocyte-derived) Gene

Reference: PTXBC041595

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FARP1-IT1; PLEKHC2; PPP1R75; FERM, RhoGEF and pleckstrin domain-containing protein 1; FERM, RhoGEF (ARHGEF) and pleckstrin domain protein 1 (chondrocyte-derived); PH domain-containing family C member 2; chondrocyte-derived ezrin-like protein; pleckstrin h		

For research use only.

GENE DESCRIPTION	FERM, RhoGEF (ARHGEF) and pleckstrin domain protein 1 (chondrocyte-derived)
GENE ID	10160
ORIGIN SPECIES	Human
ACCESSIONS	BC041595
VECTOR	pDONR223
NCBI SYMBOL	FARP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1