

C1orf210-chromosome 1 open reading frame 210 Gene

Reference: PTXBC041633

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TEMP; type III endosome membrane protein TEMP; chromosome 1 open reading frame 210		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 210
GENE ID	149466
ORIGIN SPECIES	Human
ACCESSIONS	BC041633
VECTOR	pDONR223
NCBI SYMBOL	C1orf210
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatgagacaaacaaacttggggccttcggagctccacacagcgtctgtgtgccccctggccaggcactggggctcgggcatggcctgtgctgtaggattgtgctggggctgtgtgctctctcgtcctcattgcactgtgccaatgccacctctgccgcgataccat
gccagctaccggcaccgcccactgcctgagacaggaaggggaggccgcccacagggtgctgaagatgaggatgatgatgcttcacgaggacaattacattcagcctgggactggcagctggggacagagggtagcagggaccacttccctctga