



Reference: PTXBC041801

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OFD5; ZNHIT5; DEAD (Asp-Glu-Ala-Asp) box polypeptide 59; DEAD box protein 59; zinc finger HIT domain-containing protein 5; DEAD-box helicase 59		

For research use only.

GENE DESCRIPTION	DEAD (Asp-Glu-Ala-Asp) box polypeptide 59
GENE ID	83479
ORIGIN SPECIES	Human
ACCESSIONS	BC041801
VECTOR	pDONR223
NCBI SYMBOL	DDX59
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]