

C10orf107-chromosome 10 open reading frame 107 Gene

Reference: PTXBC041932

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C10orf107; bA63A2.1; chromosome 10 open reading frame 107		

For research use only.

GENE DESCRIPTION	chromosome 10 open reading frame 107
GENE ID	219621
ORIGIN SPECIES	Human
ACCESSIONS	BC041932
VECTOR	pDONR223
NCBI SYMBOL	C10orf107
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattctctattcagttatcaaaatttatgactttactagctatgtcacttcaaaatcttaaaacactacatattgtccttagaggaaagcataaaatggcttgagagaatttatggctgaaataggaccaacacattcgcaaaagagtggagactggaatatttgatgtaaaacaagccaatg
ctatcattgattacttaaaatcagcttattcaacactacaagctatagcagtttatgttctactctccagggaagaaatgtgataggaaactgagcaagtgatagaggttgtaagctctgcatgtggccctttcccaaatcctctggaagaagggaatctattgatatttattcaacattcatagag
ccccccacaatttggatagcgaatgaagaggttgatcaagaacaagccctgaggagttctcagccagaaactgacacctcagacatggatccttttagtttggttcaccattgaagatgtaaaatcagtgctggatcaagtcacagatgacattttaatcggcattcagaccagataaac
gaaaaactgcaaatacaggaaggcccttaatgcacgaatagaaaaattgaaaaagccctaa