

C4orf32-chromosome 4 open reading frame 32 Gene

Reference: PTXBC041964

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C4orf32; chromosome 4 open reading frame 32		
For research use only.			

GENE DESCRIPTION	chromosome 4 open reading frame 32
GENE ID	132720
ORIGIN SPECIES	Human
ACCESSIONS	BC041964
VECTOR	pDONR223
NCBI SYMBOL	C4orf32
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagcactttggaagatgaagaactaaggtcttaagacgtttacatatttctcatgttatgttcattgtttgtcttttttttttttaaggaatgttgaagactcacagaaccacactggtgagccggttgagatgactacaagaaatgggaacacttttggtaactgaacaaaacctta
tcaacatgggcttcacaaggatgtattttggagaacgaatagtggaaccagtaatagtcattttctttgggttatgctgtggttcttgccctgcatgcccttgactagttgctgttctttgacctgtattattatgtgcaacagtaa