



Reference: PTXBC042054

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PUNC; HsT18880; immunoglobulin superfamily DCC subclass member 3		

For research use only.

GENE DESCRIPTION	putative neuronal cell adhesion molecule
GENE ID	9543
ORIGIN SPECIES	Human
ACCESSIONS	BC042054
VECTOR	pDONR223
NCBI SYMBOL	PUNC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr