

CCDC80-coiled-coil domain containing 80 Gene

Reference: PTXBC042105

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DRO1; SSG1; URB; okuribin; coiled-coil domain-containing protein 80; down-regulated by oncogenes protein 1; nuclear envelope protein okuribin; steroid sensitive gene 1; up-regulated in BRS-3 deficient mouse homolog; coiled-coil domain containing 80		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 80
GENE ID	151887
ORIGIN SPECIES	Human
ACCESSIONS	BC042105
VECTOR	pDONR223
NCBI SYMBOL	CCDC80
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcttctagtctggaagacggaatgtcaatcctggtatcctcccaatgtggtccatggtgattgtgtacgatttaattgattcgatgcaactcggagacagaaatggcgattcagcagtcactggggatgcgctgccagaagatgagtatgcaggctatggtaccatagttacca
ccaaggataccaggatggttaccaggatgactaccgtcatcatgagagttatcaccatggatacccttactga