

ZNF544-zinc finger protein 544 Gene

Reference: PTXBC042177

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 544		
For research use only.			

GENE DESCRIPTION	zinc finger protein 544
GENE ID	27300
ORIGIN SPECIES	Human
ACCESSIONS	BC042177
VECTOR	pDONR223
NCBI SYMBOL	ZNF544
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagcagcttctatgctgttccacccaggcatctgtgtgttcgaggatgtggtatggcattcacacaggaggagtggaacagctggacctggccagaggacactgtaccgagagtgacactggagacctgggagcatattgtctccctggggcttttcttccaaatctgat
gtgatctctcagctggagcaagaagaggacctgtgcagggcagagcaggaggccccccagggtacttcaggctcttcacaagtcacagggtgtgtggtgatgggagcttcagggtcccggtctgggacttcacataa