

ATP11B-ATPase, class VI, type 11B Gene

Reference: PTXBC042180

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	P4-ATPase flippase complex alpha subunit ATP11B; ATP1F; ATP1R; ATPase IR; ATPase, class VI, type 11B; truncated ATPase 11B protein; ATPase phospholipid transporting 11B (putative)		

For research use only.

GENE DESCRIPTION	ATPase, class VI, type 11B
GENE ID	23200
ORIGIN SPECIES	Human
ACCESSIONS	BC042180
VECTOR	pDONR223
NCBI SYMBOL	ATP11B
LENGTH	Proteogenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggcgtgatccggcagcagctgggtttgaccaccacatcagatgacacaagaaccatctacgtagccaacaggttctcagaatggcctttacacacctcagaaattatagataacagatcattctcaagtacactgtgtggaatttttccaaaaattattgaaacagtt
cagaagagtgccaacttttatttcttattatattttgttcagcttatgattgatacacctaccagtcagttaccagtggaacttcatttctttgtgataacagtaactgccataaaagcaggatagaagattgttacggcalaactcagataatgaagtaaatggagctcctgtttatgtgttc
gaagtggtggcctgtaaaaactagatcaaaaaacattcgggtgggtgataattgtcgaatagccaagatgaaattttcctgcagactgggtcctctcctcagatcgactggatgggtcctgtcagcttacaactgctagtgtggacgggaaactaacctgaaggtttgctgcatatgttg
agtattgctctgtga