

## MFSD1-major facilitator superfamily domain containing 1 Gene

Reference: PTXBC042197

|                |                                                                                                                                                              |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                          | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                         | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | SMAP4; major facilitator superfamily domain-containing protein 1; smooth muscle cell-associated protein 4; major facilitator superfamily domain containing 1 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | major facilitator superfamily domain containing 1                                                |
| <b>GENE ID</b>          | 64747                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC042197                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | MFSD1                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggaggaggaggatgaggaagcgcggcgctcctggcaggcggcctgacgagccgacagaggtgccccggcgccctgagccctgcccgcaccagtcgctggcgaccggccttttggtgctgtactgatgtcttccttggccttggcagctattttgctat  
gataatcctgctgcccctcagactcaagttaaacgagatagcaagtgaaatccacgaaattcatgctgctgtagcctggtattcttggcccaatgtagttttggtgcttttgatagaccgagttattggaatagcaggggcacaaatcatttttagctgctttgcatggacag  
gtgtttttgcccctgggtggaatattttagctttttggtgtaggaatttgaagattttagttgggattggtggcgagtccttagcagttggccagaatacatatgctgtgagctggtttaaaggcgaagaattaaacctgggtttggacttcaacttagcagctagaattggaagtacagtaaa  
catgaacctcatgggatggctgtattctaagattgaagcttttaggttctgctggtcacacaacctcgggacacactatgaltgggggtataacgtgtattcttcaataatctgtgcttggccttggccttggccttggatcagagagcagagagaatcctcataaagaacaggaaaaacag  
gtgaagtattaaataactgatgtaaaggacttctcctacccctgtggtctatattatcatctgtgctgctattatgtgtgtgttccctttattggacttgggaaagttttcttacagagaaatttggatttcttccaggcagcaagtgcattaaacagtggtgtatgtcatatcagctccatgt  
ccccgggttttggcctcctggtgataaaacagggaagaacatcatctgggttctttgcccagtagcagccactctgtgtcccatgatgctggcctttacgatgtggaaccttggattgctatgtgtcttgggacttccctactcattgcttgcctgtgcatgtgccaatgggtgcatttg  
tagttcctgaacatcagctgggaactgcataatgctcatgcagtcattcagaatcttgggttggccatcattccactgctggtatgatactgattctcgggggtattgttttgggaagtgttcttcatgacctgtgtttcttgcactttatctgtggttacttatttggatcgtgccag  
gggtgggaacctaaattattctgcaagacaaagggaagaataaaattttccatactgaatga