

C1orf87-chromosome 1 open reading frame 87 Gene

Reference: PTXBC042203

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C1orf87; CREF; carcinoma-related EF-hand protein; chromosome 1 open reading frame 87		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 87
GENE ID	127795
ORIGIN SPECIES	Human
ACCESSIONS	BC042203
VECTOR	pDONR223
NCBI SYMBOL	C1orf87
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcctgaagtccccgagatgtctcaagcaaaactgaacatatgaaactccagaagagagctgcagccagaaagctctctgtgaaacttcagcctgcaagatcctctgaaacctttaagatcagccagctctccagccctctgtgaatccagctgtgaagaacaaggctg
aggaatgtgagacgtgtagacaggttcaggaaagctggaaaatgccctctacctgtgtgatctgagtaacacaggagttcttgagaagggaacgagccagacgctcattcacaactacaatctcatttacaacctgtccctgagccctcagaaaatcgaccaggccttgccgagattccgt
tcgggagaaaatatgtcttgagccagcactgcgtacttaaggagctatga