

ATF7-activating transcription factor 7 Gene

Reference: PTXBC042363

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATFA; cyclic AMP-dependent transcription factor ATF-7; transcription factor ATF-A; activating transcription factor 7		

For research use only.

GENE DESCRIPTION	activating transcription factor 7
GENE ID	11016
ORIGIN SPECIES	Human
ACCESSIONS	BC042363
VECTOR	pDONR223
NCBI SYMBOL	ATF7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagacgacagaccgtttgtgtgcaatgccacggcgtgtggacagagattacaacgaggaccctggcagttcataaacacaagcatgagatgacattgaaattggccagcccgaaactgattcagtcattgcagatcaaacgcctactccaactagattcctgaagaactg
tgaggaggtgggactctcaatgaactagctagctcttgaacatgaattcaagaagctgcagatgaggatgagaaaaaggcaagaagcaggactgtgccaaaaactgggtgtattcagacctaggctattttattgtgcttgggataatttttaattggttaa