

## PLEKHA1-pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 1 Gene

Reference: PTXBC042458

|                |                                                                                                                                                                                                                                        |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | pleckstrin homology domain-containing family A member 1; pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 1; tandem PH domain-containing protein 1; pleckstrin homology domain containing A1 |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 1      |
| <b>GENE ID</b>          | 59338                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC042458                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | PLEKHA1                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgccttatgtggatcgcagaatgcatttgggtttctagacattgaagaaatgaaacagtgaggaaattctcgaaggtacttactatgataccagagaagatgttcgtgtgtgacatggataatccacagaacctactcttggatcatcacgtgttggagccattaagcttacac  
atttcaaggttagcgtactaagctaaagccaaagcggagttctgtttgtatgaatgcaggaaatgaggaagtacttctacaagccaatgatcagcaggacatgtggaatgggttaaatgtgttaacaaagctataaaattacagtacaaagcagtcagactcacagcctaattc  
tgataaccctaagtcgccatggtgaatgtgggaaaaagcaaggtcttacagaactgatatgtgtgtggcgtaccatcattactccactcagaagaagaagtaaatgaatgtgtgaaagtattgacagaataatctgaaacgggtcacaaagccatcttcttacttactcctaaccac  
caagatagtcggttatcaagctggtatgtgttaaaacaggagcagtgatgaaacactggaagaagatatttcaattggatgaaacacacataaggctacttcaaatctgaactgaaaggaacctcttctgtataaccactaaagaggttcataaagtcagggaatgaagcaaa  
gagacataatgatgagggaacctctttgaaattgtaacaacgtctcgaacttctatgtgcagcgtgatagccctgaagagatgcacagttggattaaagcagttcttggcgcattgtgacacagcgggtcccgccagatctgcgtctctgagcatccccccggtcctcagaatccaa  
acacgcttccgtctcctaccaacgcagccaccgcccacatccacagcctctcgcagcaactcttggctcacaactttaccatggagaagcaggaattttacagatctcttgcgaaggtcaagccagggaacttcaaggtccagactgtctcgaagagaaccagcttccaaagtga  
gaacaagctctgttaagacctcaagtaaaatggccctcaggaagaaatgtgacctagtagacttggacgatgcgagcctccggtcagtgacgtgtga