

KIF26B-kinesin family member 26B Gene

Reference: PTXBC042481

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	kinesin-like protein KIF26B; kinesin family member 26B		
For research use only.			

GENE DESCRIPTION	kinesin family member 26B
GENE ID	55083
ORIGIN SPECIES	Human
ACCESSIONS	BC042481
VECTOR	pDONR223
NCBI SYMBOL	KIF26B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgctgacgttggagcaggcagccggcagtgagcactacgacgctcgcctgctccccgccaccgctctccaacatccccacccctggtgggtcccgccagctgggtgggtccagcagccagagactgggcctttgtccccgccctgtgccactccaactacacaggt
tcgccaacaagcagcggcagcaaacccagcagccttgggtcagcaatggggcggaagaagagcgggtccccaaccaccagggccaaggtcagctccagatggccaccagtccaagcaatgggaacatcctcaattcgtggccatccaggctcaccagtacctggatggcac
ctggtccctgtcgagaaccaacgggtcacctgtaccataccaggactccgaagctcttgcacacgtggcaggggcccgttagaagactcggcggctgcgccagcgtgatgagcgaggacagggagcatga