



Reference: PTXBC042498

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAA0528; CDP138; C2 domain-containing protein 5; 138 kDa C2 domain-containing phosphoprotein; C2 calcium dependent domain containing 5		

For research use only.

GENE DESCRIPTION	KIAA0528
GENE ID	9847
ORIGIN SPECIES	Human
ACCESSIONS	BC042498
VECTOR	pDONR223
NCBI SYMBOL	KIAA0528
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]