

ZDHHC21-zinc finger, DHHC-type containing 21 Gene

Reference: PTXBC042558

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	palmitoyltransferase ZDHHC21; 9130404H11Rik; DHHC-21; DHHC21; DNZ1; HSPC097; zinc finger DHHC domain-containing protein 21; zinc finger, DHHC domain containing 21; zinc finger DHHC-type containing 21		

For research use only.

GENE DESCRIPTION	zinc finger, DHHC-type containing 21
GENE ID	340481
ORIGIN SPECIES	Human
ACCESSIONS	BC042558
VECTOR	pDONR223
NCBI SYMBOL	ZDHHC21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtccggattcactttgttgaccacatgggtgctgcatgggttgattgctttgttgggtatacaatattgtttaattcccaaatgtccttctcactatgaagaaggacatattccaggcatattaataataatattctatggcatttccatattctgtctggttgccttagtgagg
gcctccataactgatccagggaagactccctgagaaacccaagatcccatgggagaaaggagttctgggaattatgtaacaagtgttaattgatgagaccaaagcgttccatcactgcagccgctgcggccactgtgtgaggagaaatggatcactgtccatggattaacaattgtgtt
ggfagaataatcattggctctttctgcaattgtgtttctacactgaacttcttactgctacgcactgatgtttcttctgccactattactatttctccactaaaaagcgttaattggacctcttgttttagacatgaattggccataatgagactagcagcctttatgggcattactatgttagttgga
ataactggactctttacactcaactaattggcatcatcacagatacaacatctattgaaaagatgtcaaacgtgttggaagatatatcgaggccccgaaagccatggcagcagaccttccagaagtgttggcactgttgaagatcctgtgtgttattccttcaggcagaggcaaccactgc
gagttccctaccactttgccaatcatgtctaa