



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

PGC-progastricsin (pepsinogen C) Gene

Reference: PTXBC042578

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

PEPC; PGII; gastricsin; pepsin C; pepsinogen C; pepsinogen group II; preprogastricsin; progastricsin

For research use only.

GENE DESCRIPTION	progastricsin (pepsinogen C)
GENE ID	5225
ORIGIN SPECIES	Human
ACCESSIONS	BC042578
VECTOR	pDONR223
NCBI SYMBOL	PGC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagtggatgggtggcttggtcgtcctccagctctggaggcagcagtggtcaaaagtccctgaagaattaaagtatcctgagaccatgaaggagaaagggcttgctgggggagttcctgaggaccacaagaatgatcctgcttggaaagtaccgcttgggtgacctcagcgtg
acctacgagccatgctcctacatggaatgctgctcctcttggtagatcagatcgggactccaccacagaacttctgctccttttgacaccggctctccccctctcctctgacctcgacctcctag