



Reference: PTXBC042607

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PDIA7; protein disulfide-isomerase-like protein of the testis; protein disulfide isomerase family A, member 7; protein disulfide isomerase like, testis expressed		

For research use only.

GENE DESCRIPTION	protein disulfide isomerase-like protein of the testis
GENE ID	204474
ORIGIN SPECIES	Human
ACCESSIONS	BC042607
VECTOR	pDONR223
NCBI SYMBOL	PDILT
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]