



Reference: PTXBC042637

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACVRLK1; ALK-1; ALK1; HHT; HHT2; ORW2; SKR3; TSR-I; serine/threonine-protein kinase receptor R3; TGF-B superfamily receptor type I; activin A receptor type II-like 1; activin A receptor type IL; activin A receptor, type II-like kinase 1; activin A recepto		

GENE DESCRIPTION	activin A receptor type II-like 1
GENE ID	94
ORIGIN SPECIES	Human
ACCESSIONS	BC042637
VECTOR	pDONR223
NCBI SYMBOL	ACVRL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]