

SH3GL3-SH3-domain GRB2-like 3 Gene

Reference: PTXBC042864

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CNSA3; EEN-B2; HsT19371; SH3D2C; SH3P13; endophilin-A3; SH3 domain containing GRB2 like 3; SH3 domain containing GRB2 like endophilin A3; SH3 domain protein 2C; SH3 domain-containing GRB2-like protein 3; SH3-domain GRB2-like 3; endophilin-3; SH3 domain co		

For research use only.

GENE DESCRIPTION	SH3-domain GRB2-like 3
GENE ID	6457
ORIGIN SPECIES	Human
ACCESSIONS	BC042864
VECTOR	pDONR223
NCBI SYMBOL	SH3GL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgttgccggcgtgaagaagcagttccacaagccagccagctatttagtgaataaataagtggtgctgaaggaaactaaactagacgatgaattctgacatggaaaggaaaatagattaccaataaagttgtgcagaaattcttcaaaaccactgaatatctcagccaaatcc
agcatcacagagctaagctaggaatgctgaacactgtgtcgaagatccgaggcgagtgaaagaccacaggataccgcagacggaaggcttgcctgggggactgtatgctgaatacgggaaggagctcggggaagactccaccttggcaatgcatfagataagttggtgaatcatg
aagctaattggctgagtgaaagactctctgataattatgaagcaaacctttatgatccacttcagttactacaagataaagatttaaagagatcgggcatcacctgaaaaagctggaaggccgcccctggattacgattataaaagaaacagtaggtgaatataccagacgaagaag
tcagacaagcggtagaaaaattgaagagtcgaaggagttggctgaagaagcagctgttaactttttgaaaatgatgtagaacaagtcagccagttggctgtgttcatagaggcagcattagactatcacagacagtcacagagattctgcaggagctgcagagcaagctacagatgcga
atatcagctgcatccagtggtcccgacgagaataacaagccaaggcctgtgaaaaggagtttctagtgcagctcaatggagtttccaccacctgtgtagtgaagacgacagcctacagcagactggaaccagcagattag