



Reference: PTXBC042914

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SEMA-IV; SEMA4; SEMAK; semaphorin-3F; sema domain, immunoglobulin domain (Ig), short basic domain, secreted, (semaphorin) 3F; sema domain, immunoglobulin domain (Ig), short basic domain, secreted, 3F; semaphorin III/F; semaphorin IV; semaphorin 3F		

For research use only.

GENE DESCRIPTION	sema domain, immunoglobulin domain (Ig), short basic domain, secreted, (semaphorin) 3F
GENE ID	6405
ORIGIN SPECIES	Human
ACCESSIONS	BC042914
VECTOR	pDONR223
NCBI SYMBOL	SEMA3F
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1