

C21orf45-chromosome 21 open reading frame 45 Gene

Reference: PTXBC042917

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf45; B28; C21orf46; FASP1; MIS18alpha; hMis18alpha; protein Mis18-alpha; FAPP1-associated protein 1; MIS18 kinetochore protein homolog A; MIS18 kinetochore protein A		

For research use only.

GENE DESCRIPTION	chromosome 21 open reading frame 45
GENE ID	54069
ORIGIN SPECIES	Human
ACCESSIONS	BC042917
VECTOR	pDONR223
NCBI SYMBOL	C21orf45
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaggcgttcggtcactgaggtgttagcagagatgcgtggtggtgagtgccgacaaaggcgaatgcagcgactcctcgtgttgggcaagagactctccgaagactcgagccgaccagctgttcagaaagtggcgagcatgtggagctccatgagcgaagacgc
gtcgttgccgacatggagaggcgagctggaggagggcggtgctgcggagagagggcggtggttctgtgctccggtgctccggcgccgctggcgactcgtgagctgggtggccagccaggagacaccaactgcacccctgctcgtgttctctgtaatgtttc
tgttgataaggacaagctatccaaacgtgaaaagaaaatggttgcgtcttgagactttgtgctgcgcgggtgctcactcaatcttggtacgtgtacagatgcagcccaagaatcttgattacaagagagactgttttgcctcagttgtgaagccattgaaagtatgttttaggggc
ctctgaaaagcaaatgtgtcagaagataaagagcttttaattctgaaagcagagtgaaatagaaaagtcttaacacagatggaagatgtcttgaagcattacaatgaagctgtgggagggccgaatccaattgtcctttgccactgttaaaagctga