

C14orf132-chromosome 14 open reading frame 132 Gene

Reference: PTXBC042922

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C14orf132; C14orf88; chromosome 14 open reading frame 132		
For research use only.			

GENE DESCRIPTION	chromosome 14 open reading frame 132
GENE ID	56967
ORIGIN SPECIES	Human
ACCESSIONS	BC042922
VECTOR	pDONR223
NCBI SYMBOL	C14orf132
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggcgagctcttctactccaagccaacgctgtctctcccttccatgaaatcaaggtcaaggcgaataagactccctgctccactctacccccagagagaaatgattctcgtctcttcagatccccagatctgaggagaaaggatgggaggaggcagcagcattt
cgctggaaggcagcagatgctttccagccccggttcagctggaaggcttggaggccgagaccactctgctctctgaagtgggtccctggagaccgaaggctcagtgagctctgtgtgtcagcactgctgctgatccctgcaagacaaatggcacttctcttcaga
agcatcatctgccttcattattagcagtaataattatccagttattatcttaccgggtgccagttttgcacatcttttgtgctctattgtgtctcatcttctcaaatggccctgggggcagggaatgaggatgcagagagatgcacgttaa